

SARS-CoV-2 (2019-nCoV) Spike S1-His

Recombinant Protein Cat No: 91159108

Product Information:

The spike (S) glycoprotein of coronaviruses contains protrusions that will only bind to certain receptors on the host cell. Known receptors bind S1 are ACE2, angiotensin-converting enzyme 2; DPP4, dipeptidyl peptidase-4; APN, aminopeptidase N; CEACAM, carcinoembryonic antigen-related cell adhesion molecule 1; Sia, sialic acid; O-ac Sia, O-acetylated sialic acid. The spike is essential for both host specificity and viral infectivity. The term 'peplomer' is typically used to refer to a grouping of heterologous proteins on the virus surface that function together. The spike (S) glycoprotein of coronaviruses is known to be essential in the binding of the virus to the host cell at the advent of the infection process. It's been reported that SARS-CoV-2 (COVID-19 coronavirus, 2019-nCoV) can infect the human respiratory epithelial cells through interaction with the human ACE2 receptor. The spike protein is a large type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which is responsible for recognizing the cell surface receptor. S2 contains basic elements needed for the membrane fusion. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.

The main functions for the Spike protein are summarized as:

Mediate receptor binding and membrane fusion;

Defines the range of the hosts and specificity of the virus;

Main component to bind with the neutralizing antibody;

Key target for vaccine design;

Can be transmitted between different hosts through gene recombination or mutation of the receptor binding domain (RBD), leading to a higher mortality rate.

Protein Construction: A DNA sequence encoding the SARS-CoV-2 (2019-nCoV) Spike Protein (S1 Subunit) (YP_009724390.1) (Val16-Arg685) was expressed with a polyhistidine tag at the C-terminus.

Source:

2019-nCoV

Expression Host:

Baculovirus-Insect Cells

Qc Testing

Purity:

> 95 % as determined by SDS-PAGE.

Bioactivity:

Measured by its binding ability in a functional ELISA. Immobilized human ACE2 protein (Fc tag) at 10µg/mL (100µL/well) can bind 2019-nCoV Spike Protein (S1, His tag, the EC₅₀ of 2019-nCoV Spike Protein (S1, His tag) (40591-V08B1) is 300-600 ng/mL.

Endotoxin:

< 1.0 EU per µg protein as determined by the LAL method.

Predicted N Terminal:

Val16

Molecule Mass:

The recombinant SARS-CoV-2 (2019-nCoV) Spike Protein (S1 Subunit, His tag) consists of 681 amino acids and predicts a molecular mass of 76.45 kDa.

Formulation:

Lyophilized from sterile 20 mM Tris, 300 mM NaCl, 10 % glycerol, pH 7.5.

Stability & Storage:

Samples are stable for up to twelve months from date of receipt at -20°C to -80°C
Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution: Detailed reconstitution instructions are sent along with the products.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT