

Cat No: Kab00114

Product Particulars: anti-ACSS3-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): Activates acetate so that it can be used for lipid synthesis or for energy generation.

Immunogen: acyl-CoA synthetase short-chain family member 3

Synonyms: ACSS3, FLJ21963

Calculated MW: 75kDa

Uniprot ID: Q9H6R3

Specificity: Human, Mouse , Rat

Tested Application: ELISA, IHC, WB

Recommended Dilution: WB :1:200-1:1000; IHC:1:20-1:200

Gene ID: 79611

Gene Location: Mitochondrion

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT