

Cat No: Kab01674

Product Particulars: anti-CHRNA1-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): After binding acetylcholine, the AChR responds by an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane.

Immunogen: cholinergic receptor, nicotinic, alpha 1 (muscle)

Synonyms: ACHRA, CHNRA

Calculated MW:

Uniprot ID: P02708

Specificity: Human, Mouse , Rat

Tested Application: ELISA, IHC, WB

Recommended Dilution:

Gene ID: 1134

Gene Location:

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT**