

Cat No: Kab01834

Product Particulars: anti-COL4A3BP-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): Shelters ceramides and diacylglycerol lipids inside its START domain and mediates the intracellular trafficking of ceramides and diacylglycerol lipids in a non-vesicular manner.

Immunogen: collagen, type IV, alpha 3 (Goodpasture antigen) binding protein

Synonyms: CERT, STARD11

Calculated MW:

Uniprot ID: Q9Y5P4

Specificity: Human, Mouse , Rat

Tested Application: ELISA, WB, IHC

Recommended Dilution:

Gene ID: 10087

Gene Location:

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT**