

Cat No: Kab02019

Product Particulars: anti-CSN3-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): Kappa-casein stabilizes micelle formation, preventing casein precipitation in milk.

Immunogen: casein kappa

Synonyms: CASK, CSN10, CSNK

Calculated MW:

Uniprot ID: P07498

Specificity: Human

Tested Application: ELISA,IHC

Recommended Dilution:

Gene ID: 1448

Gene Location:

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT