

Cat No: Kab03439

Product Particulars: anti-GFRA2-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): Receptor for neurturin. Mediates the NRTN-induced autophosphorylation and activation of the RET receptor. Also able to mediate GDNF signaling through the RET tyrosine kinase receptor.

Isoform 2: Participates to NRTN-induced 'Ser-727' phosphorylation of STAT3.

Immunogen: GDNF family receptor alpha 2

Synonyms: GDNFRB, RETL2, TRNR2

Calculated MW:

Uniprot ID: O00451

Specificity: Human, Mouse

Tested Application: ELISA,WB,IHC

Recommended Dilution:

Gene ID: 2675

Gene Location:

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT