

Cat No: Kab00036

Product Particulars: anti-ABCD1-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): Probable transporter. The nucleotide-binding fold acts as an ATP-binding subunit with ATPase activity.

Immunogen: ATP-binding cassette, sub-family D (ALD), member 1

Synonyms: ABC42, ABCD1, Adrenoleukodystrophy protein, ALD, ALDP, AMN

Calculated MW: 80kDa

Uniprot ID: P33897

Specificity: Human

Tested Application: ELISA,WB

Recommended Dilution: WB : 1:500-1:5000

Gene ID: 215

Gene Location: Peroxisome membrane

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT