

Cat No: Kab03462

Product Particulars: anti-GIPC1-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): May be involved in G protein-linked signaling.

Immunogen: GIPC PDZ domain containing family, member 1

Synonyms: C19orf3, GIPC, RGS19IP1

Calculated MW:

Uniprot ID: O14908

Specificity: Human, Mouse , Rat

Tested Application: ELISA, WB, IHC, IP, IF

Recommended Dilution:

Gene ID: 10755

Gene Location:

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT