

Cat No: Kab03499

Product Particulars: anti-GLUD2-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): Important for recycling the chief excitatory neurotransmitter, glutamate, during neurotransmission.

Immunogen: glutamate dehydrogenase 2

Synonyms: GLUDP1

Calculated MW:

Uniprot ID: P49448

Specificity: Human, Mouse , Rat

Tested Application: ELISA, WB, IHC, IF

Recommended Dilution:

Gene ID: 2747

Gene Location:

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT