

Cat No: Kab00357

Product Particulars: anti-AMD1-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): Essential for biosynthesis of the polyamines spermidine and spermine. Promotes maintenance and self-renewal of embryonic stem cells, by maintaining spermine levels (By similarity).

Immunogen: adenosylmethionine decarboxylase 1

Synonyms: AMD1, SAMDC, AMD, ADOMETDC, FLJ26964

Calculated MW: 38kDa or 21kDa

Uniprot ID: P17707

Specificity: Human, Mouse ,Rat

Tested Application: ELISA,IHC

Recommended Dilution: IHC:1:20-1:200

Gene ID: 262

Gene Location: cytosol

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT