

Cat No: Kab00569

Product Particulars: anti-ARIH2-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): E3 ubiquitin-protein ligase mediating 'Lys-48'-and 'Lys- 63'-linked polyubiquitination and subsequent proteasomal degradation of modified proteins. May play a role in myelopoiesis.

Immunogen: ariadne homolog 2 (Drosophila)

Synonyms: ARI2, TRIAD1

Calculated MW:

Uniprot ID: O95376

Specificity: Human,Mouse ,Rat

Tested Application: ELISA,WB

Recommended Dilution:

Gene ID: 10425

Gene Location:

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT