

Cat No: Kab08851

Product Particulars: anti-TOB2-antibody

Pack Size: 100µg

Pack Size: Pack Size: 100µg

Purify: Immunogen affinity purified

Host: Rabbit

Isotype: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background (Function): Anti-proliferative protein inhibits cell cycle progression from the G0/G1 to S phases.

Immunogen: transducer of ERBB2, 2

Synonyms: KIAA1663, TOB4, TROB2

Calculated MW:

Uniprot ID: Q14106

Specificity: Human, Mouse , Rat

Tested Application: ELISA, WB, IHC, IF

Recommended Dilution:

Gene ID: 10766

Gene Location:

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**OUR KINESISDX COMMITMENT : GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT**