

Catalogue No: P8053

Size: 100 ug

Product Name: Rotavirus VP7 Antibody

Species: Rotavirus A

Raised In: Rabbit

Immunogen: Recombinant Rotavirus A VP7 protein

Clonality: Polyclonal

Isotype: IgG

Form: liquid, Preservative: 0.03% Proclin 300, 50% Glycerol, 0.01M PBS, PH 7.4

Purification: >95%, Protein G purified

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

ELISA 1:10,000 - 1:20,000

Not tested in other applications.

Storage: Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

OUR KINESISDX COMMITMENT: GUARANTEED QUALITY WITH EXPERT TECHNICAL SUPPORT